



7700/7800/7900/8800/8900 - Sampler Cones

Sampler Cones




| | | |
|---|---------------|---|
| Nickel Sampler | AT7701-Ni |  |
| Nickel Plated Sampler | AT7701-Ni/Ni |  |
| Aluminium Sampler | AT7701-Al | |
| Platinum Sampler | AT7706-Pt | |
| Platinum Sampler (18mm insert) | AT7706A-Pt |  |
| Platinum Sampler Boron free | AT7706-Pt-BF | |
| Nickel Plated Sampler | AT7706-Pt/NiP | |
| Platinum Sampler with Nickel Base | AT7706-Pt/Ni | |
| Platinum Sampler with Nickel Base (15mm insert) | AT7706B-Pt/Ni | |

7700x/7800/8800 - Skimmer Cones

Skimmer Cones

| | | |
|--|---------------|---|
| Nickel Skimmer | AT7702X-Ni |  |
| Aluminium Skimmer | AT7702X-Al | |
| Platinum Skimmer with Copper Base | AT7708X-Pt |  |
| Platinum Skimmer Cone with Copper Base | AT7708X-Pt-BF | |
| Platinum Skimmer Cone with Nickel Base | AT7708X-Pt/Ni |  |

7700s - Skimmer Cones



| | | |
|--|---------------|---|
| Nickel Skimmer | AT7702S-Ni |  |
| Platinum Skimmer Cone | AT7708S-Pt |  |
| Platinum Skimmer Cone with Nickel Base | AT7708S-Pt/Ni |  |

7900/8900 - Skimmer Cones

| | | |
|--|---------------|---|
| Nickel Skimmer Cone | AT7902X-Ni |  |
| Aluminium Skimmer Cone | AT7902X-Al | |
| Platinum Skimmer Cone with Copper Base | AT7908X-Pt |  |
| Platinum Skimmer Cone with Nickel Base | AT7908X-Pt/Ni |  |

Cone Options

ConeGuard Thread Protectors

| | | |
|---|-------------|---|
| Skimmer Cone Guard For 7700s/7900 | 70-803-1004 |  |
| Skimmer Cone Guard For 7700x/7800/8800/8900 | 70-803-1008 |  |

Retaining Ring AT7704
For 7700 Sampler Cone



Graphite Gasket (PKT 3) AT7703
For 7700/7900/8800/8900 Sampler Cone

Long Life Shield Plate AT5004
For 7700/7800/7900/8800/8900

