

## Glass DC Nebulizers

**Conikal**  
1mL/min uptake      A41-1-UC1S

**SeaSpray**  
2mL/min uptake      A41-1-USS2S

**Slurry**  
4mL/min uptake      A41-1-US6

**MicroMist**  
600µL/min uptake      A41-07-UM06S  
200µL/min uptake      A41-07-UM02  
100µL/min uptake      A41-07-UM01  
50µL/min uptake      A41-07-UM005



## DC Nebulizers for HF

**OpalMist**  
1mL/min uptake      A41-07-PFA1S  
600µL/min uptake      A41-07-PFA06S  
200µL/min uptake      A41-07-PFA02  
100µL/min uptake      A41-07-PFA01  
50µL/min uptake      A41-07-PFA005



**DuraMist**  
1mL/min uptake      A41-07-DM1S  
400µL/min uptake      A41-07-DM04



**VeeSpray**  
0.6 - 3mL/min uptake      A41-1-CV6



## IsoMist XR

### Programmable Temperature Cyclonic Spray Chamber

Glass spray chamber      KT-1067-XR  
Quartz spray chamber      KT-1067Q-XR  
PFA spray chamber      KT-1067P-XR



## Spray Chambers

**Tracey**  
50mL Cyclonic      20-808-8881HE  
With B14 gravity drain



**Twister**  
50mL Cyclonic      20-809-9119HE  
With B14 gravity drain



**Twister**  
50mL Cyclonic      20-809-0222HE



**Cinnabar**  
20mL Cyclonic      20-809-0164HE



**Drain Device**      21-808-8214



### Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch

Glass Transfer Tube 86mm      21-809-2846



## Spray Chambers for HF

### Tracey TFE

50mL Cyclonic 20-809-2508  
 With B14 gravity drain  
 Suitable for use with HF



### HF Drain Device

21-809-2540



### Mounting bracket for HF Spray Chamber

21-809-2839



### Transfer Tube, Spray Chamber to ICPE-9000/9800 Torch

PFA Transfer Tube 82.5mm 31-808-3071

## Torches

### Semi Demountable Torch (HF)

30-808-8039

Suitable for ICPS-7500/8100 and ICPE-9000/9800

#### Comprises:

Torch Body 31-808-1258  
 Torch Adaptor 31-808-1260  
 Tapered Alumina Injector 1.5mm 31-808-0013  
 Joint Clip B14-Nickel KCM-14  
 GazFit Connectors (PKT 4) for 4mm OD side arms GAZ-04



### Tapered Injectors

Alumina Injector 1.0mm 31-808-0835  
 Alumina Injector 1.5mm 31-808-0013  
 Alumina Injector 1.8mm 31-808-0015  
 Alumina Injector 2.0mm 31-808-0027  
 Alumina Injector 2.4mm 31-808-0385



## Torches

### Capillary Injectors

Sapphire Injector 2.0mm 31-808-3001

### Quartz Torches


Quartz Torch Standard 1.5mm Injector (Suitable for 7500/8100) 30-808-1179

Quartz Torch for Organics, 1.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0519

Quartz Torch for High TDS, 1.8mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0518

Quartz Torch for High TDS, 2.0mm Injector (Suitable for 7500/8100/9000/9800) 30-807-0564

Quartz Mini Torch 1.5mm Injector (Suitable for ICPE-9000/9800) 30-807-0534



## RF Coils

Copper 70-900-6000C  
 Silver 70-900-6000S  
 Gold 70-900-6000G

RF Coil Installation Kit 70-900-6000T

## Common Accessories

### Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1312



### Gaurdian In-Line Filters

Guardian in-line non-return gas filter 70-803-1362  
 Guardian in-line sample filter 70-803-1108



### Trident Internal Standard Kit

Standard Kit 60-808-1185  
 Standard Kit for HF 60-808-1150

