

Glass DC Nebulizers

SeaSpray

2mL/min uptake

A31-07-USS2

MicroMist

600µL/min uptake

A31-07-UM06

Conikal

2mL/min uptake

A31-07-UC2

Slurry

4mL/min uptake

A31-07-US6



DC Nebulizers for HF

DuraMist

1mL/min uptake

A31-07-DM1



OpalMist

2mL/min uptake

A31-07-PFA2



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

IsoMist XR Kit

KT-1015-XR

Comprises:

Glass Spray Chamber

21-809-3186

Assembly (ROL)

20-809-5176

PuraFlex Transfer Tube

21-809-3455

UniFit Right Angle Connector (PKT 3)

GAZ-04-90



Spray Chambers

Tracey

50mL Cyclonic

20-809-2665HE



Twister

50mL Cyclonic

20-809-2951HE



Adaptor

Right Angle With Ball Joint

21-809-2761

Right Angle With Ball Joint (105mm long)

21-809-4237



HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

Thermo Duo 7600 Hydride Kit

KT-1156

Thermo Duo 6000/7200/7400 Hydride Kit

KT-1179



Spray Chambers for HF

Tracey TFE

50mL Cyclonic

20-809-3469



Torches

D-Torch

30-808-2844

**Comprises:**

Base and Inner Tube

31-808-2932

Quartz Outer Tube

31-808-2840

**D-Torch Options:**

Ceramic Outer Tube

31-808-3834

Optic Fibre 170mm

31-808-2933

Optic Fibre 450mm

31-808-3117

Retaining Ring

31-808-2707

Injector Adaptor

31-808-2895

Injector Ferrule

31-808-2897



PTFE Coated

Viton O-rings (PKT 10)

VT-007

Viton O-rings (PKT 10)

V-026

O-ring kit for Thermo
6000/7000/PRO D-Torch

70-803-1921

O-ring and Seal kit for Thermo
Torch & Windows

70-803-1550

**Ceramic Outer Tube and
Optic Fibre for D-Torch**

31-808-3079

Comprises:

Ceramic Outer Tube

31-808-3834

Optic Fibre 170mm

31-808-2933

Optic Fibre 450mm

31-808-3117

**Quartz Torch**

EMT Design

31-808-2811

**Tapered Injectors for
D-Torch & EMT Torch**

Quartz Injector 1.0mm

31-808-2833

Quartz Injector 1.5mm

31-808-2834



Torches

(continued)

**Capillary Injectors for
D-Torch & EMT Torch**

Quartz Injector 2.0mm

31-808-2836

Quartz Injector 2.5mm

31-808-3246

Quartz Injector 3.0mm

31-808-3247

Alumina Injector 1.0mm

31-808-3279

Alumina Injector 2.0mm

31-808-2835

Sapphire Injector 2.0mm

31-808-2977

**Duo Periscope Window**

31-808-2677

O-ring for Periscope Window
(PKT 10)

70-803-1603

**Quartz Window**

70-803-1116

**Purged Optical Path
(POP) Assembly**

70-803-1571

Window & O-Ring for Purged
Optical Path (POP)

70-803-1892

O-ring for Thermo ICP Duo Periscope
Window, (PKT 10)

70-803-1603



RF Coils

RF Coil for 6000

Silver, Teflon Coated

70-900-4002S

Gold, Teflon Coated

70-900-4002G

**RF Coil for 7000/6000 Mark II**

Silver, Teflon Coated

70-900-4005S

Gold, Teflon Coated

70-900-4005G

**RF Coil Installation Kit**

70-900-4002T



Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1268



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108

Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes (See page 22-23 for complete list of probes)

For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler
Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon
Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon
Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon
Fibre Probe 1.0mm ID 70-803-0793



For Cetac ASX-7400, 7600 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-2940

