

Glass DC Nebulizers

MicroMist

100µL/min uptake

A31-1-UM01

50µL/min uptake

A31-1-UM005



DC Nebulizers for HF

OpalMist

100µL/min uptake

A31-1-PFA01

50µL/min uptake

A31-1-PFA005Q



Spray Chambers

High Efficiency Sample Introduction System

HE-SIS for Thermo Neoma MC-ICP-MS

KT-1172



Torches

Quartz Torch

Base for Q Torch

31-808-3494



Quartz Tube Set

31-808-3357



Ceramic Torch

Ceramic Tube Set

31-808-4048



Torch Options:

Quartz Injector 1.0mm

31-808-3440

Quartz Injector 2.0mm

31-808-3898

Quartz Injector 2.5mm

31-808-3428

Sapphire Injector 2.0mm

31-808-4014

Tapered Alumina Injector 2.0mm

31-808-4012

O-ring Kit for Q Torch base (31-808-3494)

70-803-1535



Cones

Sampler Cones

Nickel Jet Sampler

TF1005-Ni

Nickel Sampler with Copper Core

TF1001-Ni/Cu

Platinum Sampler with Copper Core

TF1006-Pt/Cu



Skimmer Cones

Nickel Skimmer

TF1002A-Ni



Nickel X Skimmer with Copper Core

TF1002X-Ni



Platinum Skimmer (H, boron-free)

TF1008-Pt



Graphite Sampler Cone Gasket

TG5001



Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier 70-803-1268



Guardian In-Line Filters

Guardian DC in-line non-return gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes (See page 22-23 for complete list of probes)

For Cetac ASX-112FR

Guardian Probe ASX-112FR 70-803-2029



For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler
Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon
Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon
Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon
Fibre Probe 1.0mm ID 70-803-0793

