

Glass DC Nebulizers

SeaSpray

400µL/min uptake A13-1-USS04

MicroMist

400µL/min uptake A13-1-UM04

200µL/min uptake A13-1-UM02

100µL/min uptake A13-1-UM01

50µL/min uptake A13-1-UM005

Conikal

1mL/min uptake A13-1-UC1



Quartz DC Nebulizers

SeaSpray

1mL/min uptake A13-1-QSS1



DC Nebulizers for HF

OpalMist

400µL/min uptake A13-1-PFA04

200µL/min uptake A13-1-PFA02

100µL/min uptake A13-1-PFA01

50µL/min uptake A13-1-PFA005



DC Nebulizers for HE-SIS

MicroMist HE

200µL/min uptake C13-1-UFT02

200µL/min uptake E13-1-UFT02



Peltier Controlled Cyclonic Spray Chambers

Glass Spray Chamber KT-1116

Quartz Spray Chamber KT-1116Q

PFA Spray Chamber KT-1116P



Spray Chambers

High Efficiency Sample Introduction System

HE-SIS for Agilent ICP-MS KT-1155



*Twister

50mL Cyclonic with Aux. 20-809-2574HE

*Requires adaptor KT-1104.



Quartz Scott Style

For 7900 UHMI 20-809-4144



Spray Chambers for HF

*Tracey PFA44

44ml Cyclonic with Aux. 20-809-2578

Spray Chamber Adaptor KT-1104

*Requires adaptor KT-1104.



Torches

Quartz Torches

with 1.0mm Injector 30-807-0558

with 1.5mm Injector 30-807-0557

with 2.0mm Injector 30-807-3608

with 2.5mm Injector 30-807-0556



Torch Options

Quartz Torch Bonnet 31-808-2577

Quartz Connector with Dilution Port 31-808-3217

Quartz Connector with 45deg Dilution Port 31-808-3280

Straight Quartz Connector 31-808-3288

Clamp (PKT 2) 70-803-0846



Torches

(continued)

Semi-Demountable Torch

with 2.4mm Quartz Injector 30-808-3853



Torch Options

Semi-Demountable Torch body

31-808-3857



1.0mm Tapered Quartz Injector

31-808-4172

1.5mm Tapered Quartz Injector

31-808-4173



2.4mm Tapered Quartz Injector

31-808-4168

1.0 mm Tapered Alumina Injector

31-808-4169



1.5 mm Tapered Alumina Injector

31-808-4170

2.4mm Tapered Alumina Injector

31-808-3855

2.4mm Sapphire Capillary Injector

31-808-4159



PEEK Injector Adaptor

31-808-4478

Viton O-rings for FDT/SDT ball joints (PKT 10)

70-V-011



RF Coils

Copper

70-900-7800C



Silver

70-900-7800S



Gold

70-900-7800G



RF Coil Installation Kit

70-900-7500T




Cones

Sampler Cones


Nickel Sampler

AT7701-Ni



Nickel Sampler with Nickel plated base

AT7701-Ni/Ni



Aluminium Sampler

AT7701-Al



Platinum Sampler

AT7706-Pt

Platinum Sampler (18 mm Insert)

AT7706A-Pt



Platinum Sampler with Nickel plated base

AT7706-Pt/NiP




Platinum Sampler with solid Nickel base

AT7706-Pt/Ni

Platinum Sampler (15mm insert) with Nickel base

AT7706B-Pt/Ni



Platinum Sampler with Nickel plated base (Extended Life)

AT7706L-Pt/NiP

Skimmer Cones

Nickel Skimmer for x-lens

AT7902X-Ni



Aluminium Skimmer for x-lens

AT7902X-Al



Platinum Skimmer for x-lens

AT7908X-Pt




Cones

(continued)

Skimmer Cones (continued)

Platinum Skimmer with
Nickel base for x-lens AT7908X-Pt/Ni



Nickel Skimmer for s-lens AT7702S-Ni

Copper Skimmer for s-lens AT7702S-Cu

Platinum Skimmer for s-lens AT7708S-Pt



Platinum Skimmer with
Nickel base for s-lens AT7708S-Pt/Ni

Platinum Skimmer with
Nickel base for m-lens AT8908M-Pt/Ni



Cone Options

Skimmer Cone Guard 70-803-1004




Long Life Pt Shield Plate AT5004

Graphite Gasket (PKT 3) AT7703

Retaining Ring AT7704

Brass Skimmer Adaptor
base for m lens AT7905M-BR



Brass Skimmer Adaptor
base for x lens AT7905X-BR

Stainless Steel Skimmer
Adaptor base for x lens AT7905X-SS



Captive Screw M4 for
Skimmer Base (PKT 3) 70-803-2061

Common Accessories

Elegra Argon Humidifier

Elegra Dual Argon Humidifier 70-803-1273



Elegra Stand 70-803-1581



Common Accessories

(continued)

Guardian In-Line Filters

Guardian DC in-line non-return
gas filter 70-803-1942

Guardian in-line sample filter 70-803-1108


Eluo Adaptor for Inline Filter 70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes (See page 20-21 for complete list of probes)

For Agilent I-AS

PTFE Encapsulated Carbon Fibre Probe
0.25mm ID 60-703-1010

PTFE Encapsulated Carbon Fibre Probe
0.3mm ID 60-703-1009



For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-1803

Polyimide sheathed Autosampler
Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon
Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon
Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon
Fibre Probe 1.0mm ID 70-803-0793



For SPS 3/SPS 4

Guardian Probe 0.75mm
Probe Connecting Line (Red) 70-803-1957

PTFE Encapsulated Carbon
Fibre Probe 0.75mm ID 70-803-0908

PTFE Encapsulated Carbon
Fibre Probe 1.0mm ID 70-803-0853

