

Glass DC Nebulizers

SeaSpray

2mL/min uptake

A21-1-USS2

MicroMist

400µL/min uptake

A21-1-UM04

Slurry

4mL/min uptake

A21-1-US6

Conikal

2mL/min uptake

A21-1-UC2



DC Nebulizers for HF

DuraMist

1mL/min uptake

A21-1-DM1

400µL/min uptake

A21-1-DM04

**VeeSpray**

0.6-3mL/min uptake

A21-1-CV6



IsoMist XR

Programmable Temperature Cyclonic Spray Chamber

IsoMist XR Kit

KT-1141-XR

Comprises:

Glass Twister Spray Chamber 21-809-3185

Assembly (LOL) 20-809-5175

Adaptor 12/8 for SpectroBlue 31-808-3489

Glass Ball Cup Adaptor S19 21-809-4101

Tygon Tubing 1/4 inch x 1m 1921-1M

Metal Ball Joint Clip JC-19

GazFit Connectors for 8mm OD arm (PKT 2) GAZ-08

UniFit Right Angle Connector (PKT 3) GAZ-04-90



Spray Chambers

Tracey

50mL Cyclonic

20-809-2665HE

**Twister**

50mL Cyclonic

20-809-2951HE

**Scott Spray Chamber**

20-809-4053

HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

Hydride Kit

KT-1160

**Quartz Adaptor**

21-809-3772

**Quartz Adaptor with Aux Port**

21-809-3771



Spray Chambers for HF

Tracey TFE

50mL Cyclonic

20-809-3717

Suitable for use with HF.

Requires Transfer Tube 21-809-3726.

**PuraFlex Transfer Tube**

21-809-3726



Torches

D-Torch Blue EOP 30-808-3600**D-Torch Blue TI** 30-808-3629**Comprises:**

Base and Inner Tube 31-808-3605

Quartz Outer Tube 31-808-3630

Injector Ferrule 31-808-2962

Retaining Ring 31-808-3595

Adaptor for D-Torch B19 31-808-3651

Adaptor S19-B19-90deg 21-809-4045

Metal Ball Joint Clip JC-19

D-Torch Options:

Sheath Gas Adaptor 21-809-4082

Aux Gas Line 70-803-1424

**Tapered Injectors for Blue EOP/TI****D-Torch**

Quartz Injector 1.0mm 31-808-0009

Quartz Injector 1.2mm 31-808-3643

Quartz Injector 1.5mm 31-808-0007

Quartz Injector 1.8mm 31-808-8146

Quartz Injector 2.0mm 31-808-0008

Quartz Injector 2.4mm 31-808-3541

Quartz Injector 2.5mm 31-808-3645

Alumina Injector 1.0mm 31-808-0835

Alumina Injector 1.5mm 31-808-0013

Alumina Injector 1.8mm 31-808-0015

Alumina Injector 2.0mm 31-808-0027

**Capillary Injectors for Blue EOP/TI****D-Torch**

Alumina Injector 1.8mm 31-808-0087



Torches

(continued)

**D-Torch
ARCOS II EOP**

30-808-3657

Comprises:

Base and Inner Tube 31-808-3660

Quartz Outer Tube 31-808-3658

Injector Ferrule 31-808-2962

Retaining Ring 31-808-3595

Adaptor for D-Torch B19 31-808-3651

Adaptor S19-B19-90deg 21-809-4045

Metal Ball Joint Clip JC-19

D-Torch Options:

Sheath Gas Adaptor 21-809-4082

Adaptor S19-B13-90deg 21-809-4048

**Tapered Injectors for ARCOS II EOP****D-Torch**

Quartz Injector 2.5mm 31-808-3676

Alumina Injector 2.0mm 31-808-3699

Alumina Injector 2.4mm 31-808-3700

**Capillary Injectors for ARCOS II EOP****D-Torch**

Sapphire Injector 2.0mm 31-808-3718

Sapphire Injector 2.4mm 31-808-3850

**Quartz Torch**

SpectroBlue EOP 31-808-3394

SpectroBlue TI 31-808-3468

**Injectors for Quartz Torch**

Quartz Injector 1.2mm 31-808-3410

Quartz Injector 1.8mm 31-808-3391

Quartz Injector 2.0mm 31-808-3411

Quartz Injector 3.0mm 31-808-3412

Alumina Injector 1.8mm 31-808-3409

Alumina Injector 2.5mm 31-808-3390



Common Accessories

Elegra Argon Humidifier

Elegra Argon Humidifier	70-803-1269
Elegra Stand	70-803-1581




Guardian In-Line Filters

Guardian DC in-line non-return gas filter	70-803-1942
Guardian in-line sample filter	70-803-1108
Eluo Adaptor for Inline Filter	70-803-1160



Trident Internal Standard Kits

Trident CT Standard Kit	60-703-1179
Standard Kit for HF	60-808-1150



Magnifier Inspection Tool 70-803-1923



Autosampler Probes (See page 20-21 for complete list of probes)

For Cetac ASX-200/500/800 Series

Guardian Probe 0.75mm	
Probe Connecting Line (Red)	70-803-1803
Polyimide sheathed Autosampler Probe 0.5mm ID	60-808-1186L
PTFE Encapsulated Carbon Fibre Probe 0.5mm ID	70-803-0784
PTFE Encapsulated Carbon Fibre Probe 0.75mm ID	70-803-1443
PTFE Encapsulated Carbon Fibre Probe 1.0mm ID	70-803-0793



For Cetac ASX-7400, 7600 Series

Guardian Probe 0.75mm	
Probe Connecting Line (Red)	70-803-2940

