

## Glass DC Nebulizers

### SeaSpray

2mL/min uptake	A31-07-USS2
400µL/min uptake	A31-07-USS04

### MicroMist

600µL/min uptake	A31-07-UM06
200µL/min uptake	A31-07-UM02
100µL/min uptake	A31-07-UM01
50µL/min uptake	A31-07-UM005



### Conikal

2mL/min uptake	A31-07-UC2
1mL/min uptake	A31-07-UC1

### Slurry

4mL/min uptake	A31-07-US6
----------------	------------

## DC Nebulizers for HF

### DuraMist

1mL/min uptake	A31-07-DM1
400µL/min uptake	A31-07-DM04



### OpalMist

2mL/min uptake	A31-07-PFA2
600µL/min uptake	A31-07-PFA06
200µL/min uptake	A31-07-PFA02
100µL/min uptake	A31-07-PFA01
50µL/min uptake	A31-07-PFA005



### VeeSpray

0.6 - 3mL/min uptake	A31-07-CV6
----------------------	------------



## IsoMist XR

### Programmable Temperature Cyclonic Spray Chamber

Glass Spray chamber	KT-1166-XR
---------------------	------------



## Spray Chambers

### HydraMist Simultaneous Cold Vapor/ Pneumatic Nebulization Spray Chamber

Hydride Kit	KT-1160
-------------	---------



## Spray Chambers

### Tracey

50mL Cyclonic	20-809-2665HE
---------------	---------------



### Twister

50mL Cyclonic	20-809-2951HE
---------------	---------------



### Straight Ball Joint Adaptor

21-809-2723
-------------



## Spray Chambers for HF

### Tracey TFE

50mL Cyclonic	20-809-9198
---------------	-------------



## Torches

### D-Torch

30-808-2877
-------------



### Comprises:

Base and Inner Tube	31-808-2932
Quartz Outer Tube	31-808-2878

### D-Torch Options:

Ceramic Outer Tube	31-808-3836
--------------------	-------------



PTFE Coated	
Viton O-rings (PKT 10)	VT-007
Viton O-rings (PKT 10)	V-026
O-ring kit	70-803-1520
Retaining Ring	31-808-2707
Injector Adaptor	31-808-2895
Injector Ferrule	31-808-2897



**Torches**

**Quartz Torch**

EMT Design 31-808-2810



**Tapered Injectors for D-Torch & EMT Torch**

Quartz Injector 0.75mm 31-808-3511  
 Quartz Injector 1.0mm 31-808-2833  
 Quartz Injector 1.5mm 31-808-2834



**Capillary Injectors for D-Torch & EMT Torch**

Quartz Injector 1.0mm 31-808-3510  
 Quartz Injector 2.0mm 31-808-2836  
 Quartz Injector 2.5mm 31-808-3246  
 Quartz Injector 3.0mm 31-808-3247  
 Alumina Injector 1.0mm 31-808-3279  
 Alumina Injector 2.0mm 31-808-2835  
 Sapphire Injector 2.0mm 31-808-2977



**Common Accessories**

**Elegra Argon Humidifier**

Elegra Argon Humidifier 70-803-1303




**Guardian In-Line Filters**

Guardian DC in-line non-return gas filter 70-803-1942  
 Guardian in-line sample filter 70-803-1108



**Trident Internal Standard Kits**

Trident CT Standard Kit 60-703-1179  
 Standard Kit for HF 60-808-1150



**Autosampler Probes**

(See page 15 for complete list of probes)

**For Cetac ASX-200/500/800 Series**

Polyimide sheathed Autosampler Probe 0.5mm ID 60-808-1186L  
 PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-0784  
 PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-1443  
 PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0793

