

ICP-MS Cones

Guaranteed performance
Immediate delivery

Glass Expansion has over 30 years of experience in manufacturing ICP-MS sample introduction components and has earned a world-wide reputation for quality and reliability.

We are proud to be the first and only ICP-MS cone manufacturer with the required technology and equipment to produce the full range of ICP-MS cones entirely in-house. Glass Expansion's manufacturing capabilities for this include: CNC machining, precision grinding, electron beam welding and electrodeless nickel plating. This results in our company being able to maintain the highest level of precision, quality and reliability standards.

We supply cones for all of the popular ICP-MS models and, to make sure you get all the support you need, we have a staff of technical experts and a fully equipped ICP-MS laboratory for testing and evaluation.

No-Risk Guarantee

As with all of our products, our cones are backed by our complete satisfaction guarantee.

Free Refurbishment

Extend the life of your platinum cones with our free refurbishment program. In many cases, a platinum cone can be refurbished multiple times and its lifetime greatly extended. As a customer service, Glass Expansion will refurbish your platinum cones free-of-charge.



Platinum Reclaim

If your platinum cone cannot be refurbished, we will provide a rebate for the value of platinum in the cone.

ConeGuard

The ConeGuard Thread Protector protects the cone thread from corrosive cleaning solutions and greatly extends cone life.

Contact us at geusa@geicp.com to check the availability of cones and ConeGuards for your ICP-MS model.



GLASS EXPANSION
Quality By Design

International

6 Central Boulevard
Port Melbourne
Vic 3207, Australia

Telephone: +61 3 9320 1111
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com

Americas

31 Jonathan Bourne Drive, Unit 7
Pocasset, MA 02559, USA

Telephone: 508 563 1800
Toll Free (US): 800 208 0097
Facsimile: 508 563 1802
Email: geusa@geicp.com